

CATALOG # : 2854

RESPIRATORY

Description:

Softech® Plus ETCO<sub>2</sub> Sampling Cannula, Pediatric, 50/CS (Individually Packaged)



United States

To order this product, please call:  
**866-246-6990**

Lat  
Am

- Details:
- Pediatric Cannula with 14 ft oxygen supply tubing and 14 ft CO<sub>2</sub> sampling line with a male luer connector

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.