


Home // Airway Management // Laryngoscopes // Reusable Laryngoscopes // Rüsch® GreenSpec™ Fiber Optic Reusable Laryngoscope Handles // 4411100

CATALOG # : 4411100

## RUSCH® GREENSPEC™ FIBER OPTIC REUSABLE LARYNGOSCOPE HANDLE

**Description:**

Rusch® GreenSpec™ Fiber Optic Reusable Laryngoscope Handle



© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.