

Home // Airway Management // Endotracheal Tubes // Specialty Intubation // Sheridan® Spiral-Flex® Murphy Uncuffed ETT // 5-12909

CATALOG # : 5-12909

## SHERIDAN® SPIRAL-FLEX® UNCUFFED ETT

### Description:

Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 4.5, 5/BX



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.