

Home // Airway Management // Endotracheal Tubes // Endobronchial Tubes // Sheridan® Sher-I-Bronch® Left-Sided Endobronchial Tube // 5-16028

CATALOG # : 5-16028

## SHERIDAN SHER-I-BRONCH LEFT-SIDED ENDOBRONCHIAL TUBE

### Description:

Sheridan Sher-I-Bronch Left-Sided Endobronchial Tube, 28 Fr



To order this product, please call:  
**+1 (919) 433-4999**

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.