

Home // Airway Management // Endotracheal Tubes // Endobronchial Tubes // Sheridan® Sher-I-Bronch® Left-Sided Endobronchial Tube // 5-16037

CATALOG # : 5-16037

SHERIDAN SHER-I-BRONCH LEFT-SIDED ENDOBRONCHIAL TUBE

Description:

Sheridan Sher-I-Bronch Left-Sided Endobronchial Tube, 37 Fr



To order this product, please call:

+1 (919) 433-4999

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.