

Home // Airway Management // Endotracheal Tubes // Standard Intubation // Sheridan® EZ-Endo® Endotracheal Tube // 5-22516

CATALOG # : 5-22516

SHERIDAN® EZ-ENDO® CUFFED ETT

Description:

Sheridan® EZ-Endo® Endotracheal Tube, Size 8.0, 10/BX



United States

To order this product, please call:

866-246-6990

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.