

Home // Airway Management // Resuscitation // Resuscitation // PEEP Valves // 5385

CATALOG # : 5385

PEEP VALVES

Description:

Disposable PEEP Valve – 30 mm, 10/CS (Individually Packaged)

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Details:

- Adjustable from 5 to 20 cm H₂O
- ISO standard 30 mm ID connector
- Attaches to any resuscitation bag flow diverter or filter with a 30 mm connection

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.