


CATALOG # : 5387

LIFESAVER® DISPOSABLE MANUAL RESUSCITATOR

Description:

Lifesaver® Disposable Manual Resuscitator, Adult, 6/CS (Individually Packaged)



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Details:

- Foldable resuscitation bag incorporates the same great features as the standard Lifesaver® adult bag and collapses for compact storage.
- Swivel flow diverter allows connection of PEEP valve
- Includes adult face mask (Item Code 1276)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.