


CATALOG # : 5404

HYPERINFLATION BAG SYSTEM

Description:

Hudson RCI® Hyperinflation Bag System, 1.0 L, 20/CS (Individually Packaged)



United States

To order this product, please call:

866-246-6990

Lat

Am

Details:

- Includes 7 ft Star Lumen® Oxygen Supply Tubing and 7 ft pressure tubing

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.