

Home // Airway Management // Resuscitation // Resuscitation // Hyperinflation Bag System // 5408

CATALOG # : 5408

HYPERINFLATION BAG SYSTEM

Description:

Hudson RCI® Hyperinflation Bag System, 1.0 L (Not made with natural rubber latex), 10/CS (Individually Packaged)



United States

To order this product, please call:

866-246-6990

Lat

Am

Details:

- Includes Infant Mask (Item Code 15802);
- Manometer;
- 7 ft Star Lumen® Oxygen Supply Tubing

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.