


Home // Airway Management // Laryngoscopes // Single-Use Laryngoscopes // Rüsch® DispoLED® Single-Use Fiber Optic Laryngoscope Handle // 77700

CATALOG # : 77700

RUSCH® DISPOLED® SINGLE-USE FIBER OPTIC LARYNGOSCOPE HANDLE



NO IMAGE
AVAILABLE

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.