

Home // Airway Management // Endotracheal Tubes // Endotracheal Tube Holders // Hudson RCI® Gentle-Lock™ Oral Endotracheal Tube Holder // 81801

CATALOG # : 81801

HUDSON RCI® GENTLE-LOCK™ ORAL ENDOTRACHEAL TUBE HOLDER

Description:

Endotracheal Tube Holder Standard



United States

To order this product, please call:

866-246-6990

Details:

- Standard size
- Does not include bite block
- 12/case

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.