


Catalog # :

353195

Debakey Ring Handle Bulldog Clamp

Description:

353195 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey®-type Jaws and spring handles provide firm, non-traumatizing grip. "S"-shaped, curved, 2-1/2" (6.5 cm) jaws, long and curved shafts, 5-1/8" (13 cm).



To order this product, please call:

+1 (919) 433-4999