

Catalog # :


~~507932~~**Maloney Tapered Tungsten-Filled Esophageal Bougie****Description:**

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area, cardioesophageal junction or stenosis. Once the tip is through the involved area, the tapered end allows the single bougie to dilate this portion of the esophagus to its maximum limit.

Overall length is approximately 29-1/2" (75 cm). 32 Fr.

Related:

| | | |
|---------------|---------------|---------------|
| 507920 20 Fr. | 507934 34 Fr. | 507948 48 Fr. |
| 507922 22 Fr. | 507936 36 Fr. | 507950 50 Fr. |
| 507924 24 Fr. | 507938 38 Fr. | 507952 52 Fr. |
| 507926 26 Fr. | 507940 40 Fr. | 507954 54 Fr. |
| 507928 28 Fr. | 507942 42 Fr. | 507956 56 Fr. |
| 507930 30 Fr. | 507944 44 Fr. | 507958 58 Fr. |
| 507932 32 Fr. | 507946 46 Fr. | 507960 60 Fr. |



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699