

**Catalog # :**

**518035**

**Jackson Improved Tracheostomy Tube #6**

**Description:**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator.

518035 (size 6, 10.0 mm O.D., 7.2 mm I.D., 69.0 mm long).

*Related:*

518010: Size 1

518015: Size 2

518020: Size 3

518025: Size 4

518030: Size 5

518035: Size 6

518040: Size 7

518045: Size 8

518050: Size 9

518055: Size 10

518029: Size 4 (fen.)

518034: Size 5 (fen.)

518039: Size 6 (fen.)

518044: Size 7 (fen.)

518049: Size 8 (fen.)



United States

Details:

- Stainless steel.

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator.