

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #7 FOR SHIKANI-FRENCH SPEAKING VALVE

Catalog # :


~~518040A~~

---

## JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #7 FOR SHIKANI-FRENCH SPEAKING VALVE

Description:

8.1mm I.D., 11.0mm O.D., 69mm



NO IMAGE  
AVAILABLE

---

United States

To order this product, please call:

**866-246-6990**

---

Canada

To order this product, please call:

**800-387-9699**

---