

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Esophagoscopes // ENT
// JACKSON OVAL ESOPHAGOSCOPE

Catalog # :

~~523608~~

JACKSON OVAL ESOPHAGOSCOPE

Description:

10 x 14 x 40, Overall length 401mm, min. I.D. 9.2mm, replacement light carrier 526274

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699