


Catalog # :

~~162410~~**FERGUSON FRAZIER SUCTION TUBE****Description:**

Stainless steel tubes, block-type finger valve, complete with cleaning stylette. Angled. 4-7/8" (12.5 cm). Pilling® Pattern, satin finish.

0, 6 Fr. tip

Stainless steel tubes, block-type finger valve, complete with cleaning stylette. Angled. 4-7/8" (12.5 cm). Pilling® Pattern, satin finish.



United States

Details:

- # 0, 6 Fr. tip

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699