

Catalog # :


~~166002~~**Pilling® Poly-Clamp Rail Clamp - Set of Four****Description:**

The Pilling® Poly-Clamp is a highly versatile operating table rail clamp for use with the Strong T-Bar, Gomez, Poly-Tract Fowler and other retractor systems, as well as with other equipment to be attached to the operating room table. The Poly-Clamp has been designed to fit the side rail of any currently manufactured U.S. or European operating room table.

Four clamps, each accept 5/8" (1.6 cm) round bar or flat bar uprights.

*Related:*

341114 One 341114 Pilling® Poly-Clamp Rail Clamp



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**