

**Catalog # :****341155****Chaux Sternal IMA-Mitral Valve Retractor Set****Description:**

Three retractors in one, can be used as a sternal/IMA and valve retractor. Includes four sternal blades with individual and independent dual movement. Sliding sternal blade assemblies allow movement on horizontal planes of side arms to provide proper position adjustment according to length of sternum. Both side arms can be rotated and locked at 30 degrees or 60 degrees along its axis for flexibility in positioning the retractor for IMA exposure. All components of the retractor remain in the sterile field. The valve retractor fits directly onto arms of the sternal blade assembly and provides self-retaining valve retraction. The Chaux Sternal IMA Valve Retractor when used for valve surgery consists of the basic Chaux Retractor with the addition of a Lemole-Pilling® Valve Retractor Set (341501) and Valve Sternal Blade Assembly (341149).

Spread 9-1/8" (23 cm), Side Blade 1-1/8" wide x 1-1/4" deep (2.5 x 3 cm)

**Complete set includes:**

- 341149 Valve Sternal Blade Assembly, qty. 1
- 341150 Chaux Multipurpose Sternal IMA Retractor, qty. 1
- 341511 Lemole-Pilling® T-Frame, qty. 1
- 341520 Universal Clamp for Lemole, qty. 2
- 341522 Cooley Atrial Retractor, small, 9-1/2" (24 cm), qty. 2
- 341523 Cooley Atrial Retractor, medium, 9-1/2" (24 cm), qty. 1
- 341524 Cooley Atrial Retractor, large, 9-1/2" (24 cm), qty. 1
- 341525 Malleable Atrial Retractor, 8-3/4" (22 cm), qty. 2



To order this product, please call:

**+1 (919) 433-4999**