

**Catalog # :****341200****Gerbode-Burford Retractor, Child****Description:**

Stainless steel. Has same projection-free construction as adult models but is proportionally reduced to child size.

Spread 5-15/16" (15 cm), supplied with four pairs of side blades:

- 3/4" deep x 1" wide (2 x 2.5 cm)
- 1-1/16" deep x 1-1/4" wide (2.5 x 3 cm)
- 1-1/16" deep x 1-3/4" wide (2.5 x 4.5 cm)
- 1-7/8" deep x 1-3/4" wide (5 x 4.5 cm)

Arm length 4-1/2" (11.5 cm)

Stainless steel



To order this product, please call:

**+1 (919) 433-4999**