


Catalog # :

~~341310~~**Heitmiller Mini Thoracotomy Retractor****Description:**

Designed by Dr. Richard Heitmiller.

341310 This retractor provides excellent exposure for video-assisted thoracoscopy, open lung biopsies, transaxillary incisions and mid-CAB procedures in which mini-thoracotomy incisions are used. With its small "footprint," it also is ideal for mobilizing lung adhesions during video thoracotomies.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699