

Catalog # :

~~341514~~

Chaux-Lemole Bi-Angle T-Frame for Cooley and Chaux Retractors

Description:

Bi-Angle T-Frame, for additional third point atrial retraction during mitral valve surgery, fits Cooley-style and Chaux Retractors.

Arm tip maximum diameter 0.252", overall width 9-7/8" (25 cm)

Related:

341150 Chaux Multipurpose Sternal IMA Retractor

341155 Chaux Sternal/IMA Valve Retractor Set

341194, 341195 Cooley Adult Sternum Retractors

341198 Bonchek-Carpenter Sternal/Atrial Valve Retractor

341506 Cooley Sternum Retractors



To order this product, please call:

+1 (919) 433-4999
