

Catalog # :

~~341635~~**Brewster Phrenic Retractor Nested Pair****Description:**

341635 Nested set of 2 Retractors:

Outside Retractor (double-ended) - overall length 149 mm (5-7/8")


-Large End - 13 mm wide x 35 mm deep

-Small End - 10 mm wide x 24 mm deep

Inside Retractor (double-ended) - overall length 145 mm (5-3/4")

-Large End - 13 mm wide x 32 mm deep

-Small End - 10 mm wide x 20 mm deep

A large rectangular box with a thin border, containing the text "NO IMAGE AVAILABLE" in a light gray, sans-serif font, centered within the box.

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**