

Catalog # :

~~342006~~**STERNAL LOCKING PLATES****Description:**

A device used to fixate size 5 (B&S Gage 20) monofilament stainless steel suture wire to achieve sternal approximation. Only 316L stainless steel sternal closure wire should be used with the locking plate. Suture wire of unknown origin or type may not be metallurgically compatible with the locking plates.

8 plates per package (10 packages per box). Suture not included.

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699