

Home // Cannulas, Suction Tubes, Aspirating Tubes, Needles // Suction, Aspirating Tubes, Needles // CVT // Debakey®
Heparin Needle Small

Catalog # :

~~351575~~

Debakey® Heparin Needle Small

Description:

Collectively, this cannula and those listed in "Related" below constitute a complete line of stainless steel, reusable heparin cannulae. A selection of sizes, shapes and lengths provides for a variety of uses in flushing, irrigation and injection. Applications include peripheral vascular, cardiac, general, orthopedic, urology and neuro surgery.

Debakey Heparin Needle Small, malleable, 3 mm tip, 2-1/2" (6.5 cm)

Related:

351560 Stoney Heparin Injector, 1-7/8" (4.5 cm)

351573 Zanger-Debakey Heparin Needle, malleable, 4 mm tip, 2-1/2" (6.5 cm)

351574 Horsley Heparin Needle, 2-3/8" (6 cm)

351576 Debakey Heparin Needle Medium, malleable, 3.7 mm tip, 2-7/8" (7.5 cm)

351577 Debakey Heparin Needle Large, malleable, 4 mm tip, 6" (15 cm)

351578 IMA Heparin Needle, 23 gauge tip, 2-1/8" (5.5 cm)

351586 Garrett Heparin Injector, malleable, 5/16" (8 mm) tip, working length 3" (7.5 cm)

351587 Garrett Heparin Injector, malleable, 1/2" (12 mm) tip, working length 3" (7.5 cm)

351593 Samson Heparin Injector, malleable, small 7 mm tip, 3-3/8" (8.5 cm)

351594 Samson Heparin Injector Medium, malleable, medium 9 mm tip, 3-3/8" (8.5 cm)

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699