

## Catalog # :

~~951576~~**DeBakey® Heparin Needle Medium****Description:**

Collectively, this cannula and those listed in "Related" below constitute a complete line of stainless steel, reusable heparin cannulae. A selection of sizes, shapes and lengths provides for a variety of uses in flushing, irrigation and injection. Applications include peripheral vascular, cardiac, general, orthopedic, urology and neuro surgery.

DeBakey Heparin Needle Medium, malleable, 3.7 mm tip, 2-7/8" (7.5 cm)

*Related:*

351560 Stoney Heparin Injector, 1-7/8" (4.5 cm)

351573 Zanger-DeBakey Heparin Needle, malleable, 4 mm tip, 2-1/2" (6.5 cm)

351574 Horsley Heparin Needle, 2-3/8" (6 cm)

351575 DeBakey Heparin Needle Small, malleable, 3 mm tip, 2-1/2" (6.5 cm)

351577 DeBakey Heparin Needle Large, malleable, 4 mm tip, 6" (15 cm)

351578 IMA Heparin Needle, 23 gauge tip, 2-1/8" (5.5 cm)

351586 Garrett Heparin Injector, malleable, 5/16" (8 mm) tip, working length 3" (7.5 cm)

351587 Garrett Heparin Injector, malleable, 1/2" (12 mm) tip, working length 3" (7.5 cm)

351593 Samson Heparin Injector, malleable, small 7 mm tip, 3-3/8" (8.5 cm)

351594 Samson Heparin Injector Medium, malleable, medium 9 mm tip, 3-3/8" (8.5 cm)

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**