

Catalog # :

~~351675~~**Pilling® Multipurpose Retractor Complete****Description:**

A versatile self-retaining Henly-type Retractor designed for use in the following surgical specialties:

- Peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions
- Spinal surgery in the cervical and lumbar regions
- Orthopedic surgery in shoulder and hip procedures
- General surgery for inguinal hernia, appendectomy and many minor procedures.

Ball-snap attachment of side blades. Swivel side blades, rotate 360 degrees.

Pilling® Multipurpose Retractor Set, complete with:


- 1 Frame with a 3-1/8" (7.8 cm) spread
- 3 Center Blades 7/8" (2.2 cm), 1-1/8" (2.8 cm), and 1-1/4" (3.2 cm) deep x 5/8" (1.6 cm) wide
- 3 pairs of Side Blades 1-1/8" (2.9 cm), 2-1/8" (5.5 cm), and 3-1/8" (8.0 cm) deep x 7/8" (2.2 cm) wide

Related:

Side blades for 351675 - 1-1/8", 2-1/8", 3-1/8"

Center blades for 351675 - 7/8", 1-1/8", 1-1/4", 2", 2-1/2"

Extra-long center blades for 351675



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699