

Catalog # :

~~951779~~**Henly Popliteal Retractor****Description:**

4 X 3-pronged, Supplied w/3 center blades: 5/8" (1.5cm) width x 3/4" (1.8cm), 5/8" (1.5cm) x 1" (2.5cm) and 5/8" (1.5cm) x 1 1/4" (3.1cm) depth blade, 6 3/8" (16.0cm)

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699