



Catalog # :

~~353189~~**Debakey Ring Handle Bulldog Clamp****Description:**

353189 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey®-type Jaws and spring handles provide firm, non-traumatizing grip. Angled 45 degrees, jaws 1-3/16" (3 cm) past angle, curved shafts, 5-3/8" (13.5 cm).



To order this product, please call:

+1 (919) 433-4999