


**Catalog # :****353192****Debakey Ring Handle Bulldog Clamp****Description:**

353192 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey®-type Jaws and spring handles provide firm, non-traumatizing grip. Angled 90 degrees, jaws 1-3/4" (4.5 cm) past angle, straight shafts, 4-1/8" (10.5 cm).



To order this product, please call:

**+1 (919) 433-4999**