

## Catalog # :

~~355083~~**Garrett Vascular Dilator****Description:**

Short 2.5 mm, malleable shaft

*Related:*

355080 Short 1.0 mm      355091 Long 1.0 mm

355081 Short 1.5 mm      355092 Long 1.5 mm

355082 Short 2.0 mm      355093 Long 2.0 mm

355083 Short 2.5 mm      355094 Long 2.5 mm

355084 Short 3.0 mm      355095 Long 3.0 mm

355085 Short 3.5 mm      355096 Long 3.5 mm

355086 Short 4.0 mm      355097 Long 4.0 mm

355087 Short 4.5 mm      355098 Long 4.5 mm

355088 Short 5.0 mm      355099 Long 5.0 mm

355089 Set of 9 short (1-5 mm), 5-1/2"

355090 Set of 9 long (1-5 mm), 8-1/4"

355100 Dilator holder, 8-1/8" (20.5 cm)

355101 Dilator holder, 5-1/8" (13 cm)

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**