

Catalog # :
355088

Garrett Vascular Dilator

Description:

Short 5.0 mm, malleable shaft

Related:

- | | |
|---------------------|--------------------|
| 355080 Short 1.0 mm | 355091 Long 1.0 mm |
| 355081 Short 1.5 mm | 355092 Long 1.5 mm |
| 355082 Short 2.0 mm | 355093 Long 2.0 mm |
| 355083 Short 2.5 mm | 355094 Long 2.5 mm |
| 355084 Short 3.0 mm | 355095 Long 3.0 mm |
| 355085 Short 3.5 mm | 355096 Long 3.5 mm |
| 355086 Short 4.0 mm | 355097 Long 4.0 mm |
| 355087 Short 4.5 mm | 355098 Long 4.5 mm |
| 355088 Short 5.0 mm | 355099 Long 5.0 mm |

- 355089 Set of 9 short (1-5 mm), 5-1/2"
- 355090 Set of 9 long (1-5 mm), 8-1/4"
- 355100 Dilator holder, 8-1/8" (20.5 cm)
- 355101 Dilator holder, 5-1/8" (13 cm)

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699