

Catalog # :**355500****BERRY ANGIO ACCESS TUNNELER****Description:**

Designed by William B. Berry, M.D. of Pasadena, Texas, Angio Access Tunneler (355500) is used to create forearm loop tunnels for placement of arteriovenous grafts for hemodialysis. Removable interchangeable handle and bullet tip allow:

- Simplified surgical procedure with better surgical control.
- Reduced likelihood of "kinking."
- Less patient trauma.

Diamond knurled handle for ease of control. Suture hole in bullet tip allows secure graft attachment.

Available in 10-3/4" (27.5 cm) length (105 degree arc)



To order this product, please call:

+1 (919) 433-4999