

Catalog # :

502245**Laryngoscope Holder****Description:**


Pilling® Laryngoscope Holders enable the surgeon to have both hands free during Laryngoscopy, Micro-Laryngoscopy and Laser Laryngoscopy procedures. These holders allow the surgeon to perform intricate procedures requiring both hands available for necessary instrumentation insertion through the laryngoscopes.

502245: Laryngoscope Holder , Improved Model. The external chest bar and rod easily attach to or detach from the handle of most laryngoscopes. Manually operated gear supplies leverage to the laryngoscope for exposure of the anterior commissure and allows the operator to have both hands free for instrumentation.

Related:

502265: Dedo Table-Mounted Mayo Tray

502248: Tobey Laryngoscope Holder



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699