

Catalog # :

**502248****Tobey Laryngoscope Holder****Description:**


Pilling® Laryngoscope Holders enable the surgeon to have both hands free during Laryngoscopy, Micro-Laryngoscopy and Laser Laryngoscopy procedures. These holders allow the surgeon to perform intricate procedures requiring both hands available for necessary instrumentation insertion through the laryngoscopes.

502248: Tobey Laryngoscope Holder is the newest (and most popular) version of Pilling® Laryngoscope holders. It has a quick attach/release mechanism that provides a secure hold of your laryngoscope. The Tobey eliminates the twisting and turning necessary to lock the laryngoscope into the holder. It allows for fine adjustments using the size adjustment knob for different laryngoscope handles.

**Related:**

502265: Dedo Table-Mounted Mayo Tray

502245: Laryngoscope Holder, Improved Model



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**