

Catalog # :**502260****Boston University Suspension System****Description:**

This system is used to suspend the laryngoscope during micro and laser surgery of the larynx. *Laryngoscope and Strong Laryngoscope Handle Extension (502261) are not included.*

- Supports the laryngoscope for two-handed surgery.
- Is not affected by chest movement.
- Allows use of the Boyce-Jackson position.
- Allows repositioning of the table as necessary without readjustment of the suspension system.

502260: Boston University Suspension System includes Pilling® Table Clamp (163100) and Pivoting Extension Arm Assembly (502262)

Boston University Suspension System Scope Service to adapt handle for use with BU System. All laryngoscopes must be tapped for use with the BU System. Item 502261 Strong Laryngoscope Handle Extension must be ordered along with scope for use with BU Suspension System. Contact Teleflex® Customer Service to have your laryngoscope tapped.

Related:

502265: Dedo Mayo Tray

502261: Strong Laryngoscope Handle Extension

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699