

Catalog # :

502566**Sanders Venturi Adapter For Pilling® Bronchoscopes****Description:**

Flow and pressure characteristics vary greatly according to the design of the bronchoscope or laryngoscope and any new combination should be studied carefully before clinical use.

Note: The Sanders Venturi System can be used only with tapered scopes.

Sanders Venturi Adapter, 15 mm sidearm

Related:

502562 Albert-Sanders Venturi Adapter

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699