

Home // Forceps, Hemostats, Clamps // Forceps, Hemostats // ENT // PILLING FLEXIBLE ENDOSCOPIC OVAL CUP
FORCEPS

Catalog # :

~~505674~~

PILLING FLEXIBLE ENDOSCOPIC OVAL CUP FORCEPS

Description:

The Pilling® Flexible Endoscopic Oval Cup Forceps were designed for use with Pilling® Laser and Standard Bronchoscopes. They easily pass through the instrument/fiber channel for the purpose of biopsy, foreign body retrieval and removal of char from laser ablation.

2 mm oval cup biopsy, 14-3/8" (36.5 cm)

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699