

Catalog # :


~~505690~~

PILLING FLEXIBLE ENDOSCOPIC CIRCULAR CUP FORCEPS

Description:

The Pilling® Flexible Endoscopic Circular Cup Forceps were designed for use with Pilling® Laser and Standard Bronchoscopes. They easily pass through the instrument/fiber channel for the purpose of biopsy, foreign body retrieval and removal of char from laser ablation.

2 mm circular cup biopsy, 19-5/8" (50 cm)



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699