

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Fiberoptic Cables // ENT
// Wolf to Pilling® Fiberoptic Cable Adapter

Catalog # :

50615AC

Wolf to Pilling® Fiberoptic Cable Adapter

Description:

This adapter is to be used with standard Pilling® Fiber Optic Cables P358PRA or P356PRA. Converts a Pilling® illuminator end to Wolf.

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699