

Catalog # :**508002****Polysulfone Telescope Tray**

This item has been discontinued and is no longer available.

Description:

Pilling® Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum protection and flexibility. The transparent lid will help you easily identify instruments. Includes silicone stabilizer bars to secure telescope during storage.

9-1/2" x 5-1/2" x 2-1/2" (24 x 14 x 6.5 cm)


Related:

508004 Video Laryngoscope Tray with Mat

521185 Polysulfone Storage Tray, 11" x 7"

521186 Polysulfone Storage Tray, 10-3/4" x 12-3/4"

521187 Polysulfone Storage Tray, 10-1/4" x 19-3/4"



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699**Details:**

- Pilling® Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum

protection and flexibility. The transparent lid will help you easily identify instruments. Includes silicone stabilizer bars to secure telescope during storage.

- 9-1/2" x 5-1/2" x 2-1/2" (24 x 14 x 6.5 cm)