

Catalog # :

~~508004~~**Video Laryngoscope Tray with Mat****Description:**

Pilling® Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum protection and flexibility. The transparent lid will help you easily identify instruments. Includes silicone finger mat. Tray size will accommodate laryngoscope and fiber-optic cable.

12-1/2" x 8-3/4" x 3-1/4" (31.5 x 22 x 8 cm)


Related:

508002 Polysulfone Telescope Tray

521185 Polysulfone Storage Tray, 11" x 7"

521186 Polysulfone Storage Tray, 10-3/4" x 12-3/4"

521187 Polysulfone Storage Tray, 10-1/4" x 19-3/4"



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699**Details:**

- Pilling® Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum

protection and flexibility. The transparent lid will help you easily identify instruments. Includes silicone finger mat. Tray size will accommodate laryngoscope and fiber-optic cable.

- 12-1/2" x 8-3/4" x 3-1/4" (31.5 x 22 x 8 cm)