

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #5, INNER AND OUTER TUBES AND OBTURATOR


Catalog # :

511130P

JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #5, INNER AND OUTER TUBES AND OBTURATOR

Description:

6.1mm I.D., 9.0mm O.D., 68mm, sterling silver



NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
