

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #6, INNER AND OUTER TUBES AND OBTURATOR

Catalog # :

511135P

JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #6, INNER AND OUTER TUBES AND OBTURATOR

Description:

7.1mm I.D., 10.0mm O.D., 69mm, sterling silver



To order this product, please call:

+1 (919) 433-4999
