

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #8, INNER AND OUTER TUBES AND OBTURATOR

Catalog # :

511145P

JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #8, INNER AND OUTER TUBES AND OBTURATOR

Description:

9.1mm I.D., 12mm O.D., 69mm, sterling silver

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699