

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #4 FOR SHIKANI-FRENCH SPEAKING VALVE

Catalog # :

~~518029A~~

JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #4 FOR SHIKANI-FRENCH SPEAKING VALVE

Description:

5.5mm I.D., 8mm O.D., 62mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699