

Catalog # :

~~518030~~**Jackson Improved Tracheostomy Tube #5****Description:**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator.

518030 (size 5, 9.0 mm O.D., 6.2 mm I.D., 68.0 mm long).

Related:

518010: Size 1

518015: Size 2

518020: Size 3

518025: Size 4

518029: Size 4 (fen.)

518030: Size 5

518034: Size 5 (fen.)

518035: Size 6

518039: Size 6 (fen.)

518040: Size 7

518044: Size 7 (fen.)

518045: Size 8

518049: Size 8 (fen.)

518050: Size 9

518055: Size 10

United States

Details:

- Stainless steel.

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

NO IMAGE
AVAILABLE

Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator.