

Catalog # :
~~518032~~

Non-Fenestrated Jackson Improved Inner Tube Only with 2 Shikani Valves, Size 5

Description:

Patients using standard Pilling® Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach the snap-on Shikani Valve. Each inner tube is provided with two speaking valves.

If you are using standard Jackson trachea tube 518030, then
ORDER: Shikani-French Speaking Valve 518032 (size 5, 6.2 mm I.D., 71.0 mm long).

- Related:
- 518027: Shikani-French Speaking Valve, Size 4
 - 518032: Shikani-French Speaking Valve, Size 5
 - 518037: Shikani-French Speaking Valve, Size 6
 - 518042: Shikani-French Speaking Valve, Size 7

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Details:

- Patients using standard Pilling® Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the

inner trachea tube with one of the tubes listed below, patients can attach the snap-on Shikani Valve. Each inner tube is provided with two speaking valves.