

Catalog # :

~~518034~~

Jackson Improved Tracheostomy Tube - Fenestrated #5

Description:

Patients using standard Pilling® Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach the snap-on Shikani Valve. This Jackson Tube is fenestrated to facilitate speech. The fenestration is located at a point measured from the neckplate equal to 1/3 of the length of the tube on both the inner and outer tubes. Each inner tube is provided with two speaking valves.

518034 (size 5, 9.0 mm O.D., 6.2 mm I.D., 69.0 mm long)

Related:

518010: Size 1

518015: Size 2

518020: Size 3

518025: Size 4

518030: Size 5

518035: Size 6

518040: Size 7

518045: Size 8

518050: Size 9

518055: Size 10


518029: Size 4 (fen.)

518034: Size 5 (fen.)

518039: Size 6 (fen.)

518044: Size 7 (fen.)

518049: Size 8 (fen.)



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Details:

- Patients using standard Pilling® Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the

inner trachea tube with one of the tubes listed below, patients can attach the snap-on Shikani Valve. *This Jackson Tube is fenestrated to facilitate speech. The fenestration is located at a point measured from the neckplate equal to 1/3*

of the length of the tube on both the inner and outer tubes.

- Each inner tube is provided with two speaking valves.